napari-{{cookiecutter.plugin_name}} Documentation

Release 0.1.8

{{cookiecutter.full_name}}

CONTENTS:

1	Installation	
	1.1	Stable release
	1.2	From sources
2	Indic	es and tables

CHAPTER

ONE

INSTALLATION

1.1 Stable release

To install napari-imsmicrolink, run this from terminal:

\$ pip install napari-imsmicrolink

This is the one method to install napari-imsmicrolink. You can also install from directly within napari by navigating to Plugins -> Install/Uninstall Plugins... in the menu bar.

If you don't have pip installed, this Python installation guide can guide you through the process.

1.2 From sources

The sources for napari-imsmicrolink can be downloaded from the Github repo.

You can either clone the public repository:

\$ git clone git://github.com/nhpatterson/napari-imsmicrolink

Or download the tarball:

\$ curl -OJL https://github.com/nhpatterson/napari-imsmicrolink/tarball/master

Once you have a copy of the source, you can install it with:

\$ python setup.py install

CHAPTER

TWO

INDICES AND TABLES

- genindex
- modindex
- search